


Progress Report

Kuei-Yueh Ko

01.29.2018

Experimental Process of 10X Genomics



Input

Start with hundreds to thousands of cells

Chromium System

Automated barcoding & library construction

Sequencer



Use existing short-read sequencer


10x Software Tools

Turn-key analysis & visualization

Gel bead in Emulsion droplet (GEM)

Encapsulates each tiny micro-reaction within the Chromium System





Analyzing Process of 10X Genomics

Single-Library Analysis with cellranger count

- Run **cellranger mkfastq** on the **Illumina BCL output folder** to generate **FASTQ files**.
- Run **cellranger count** on each library that was demultiplexed by cellranger mkfastq
- Optionally, run **cellranger aggr** to aggregate multiple libraries from a single experiment that were analyzed by cellranger count
- Optionally run **cellranger reanalyze** to re-run the secondary analysis on a library or aggregated set of libraries (i.e., PCA, t-SNE, and clustering).

mkfastq → **count** → **aggr** → **reanalyze**

Illumina

**Illumina
bcl2fastq**

**Binary base call
format (BCL)**

**FASTQ
File**

Bowtie / TopHat

**SAM/BAM
File**

Cufflinks

Cuffdiff

**Bowtie / TopHat / Cufflinks / Cuffdiff
Pipeline**

Bowtie/Tophat/Cufflinks/Cuffdiff RNA-Seq Pipeline

